





Microscope Objectives


Plan Achromatic
DW-PD Model



Plan Achromatic DW-PD Model			
Magnification	4X	10X	20X
Numerical Aperture	0.1	0.25	0.4
Working Distance (mm)	17.912	2.043	0.5
Slide CT (mm)	0.17	0.17	0.17
Focus Length (mm)	28.9	15.37	8.4
BFL (mm)			
Conjugate Distance (mm)	195	195	195
Tube Length (mm)	160	160	160
Wavelength (nm)	400-700	400-700	400-700
Resolution (μm)	2.8	1.1	0.7
Physical Depth of Focus (μm)	27.5	4.4	1.7
Thread	W4/5' x 1/36'		

 : www.shanghai-optics.com

 : rfq@shanghai-optics.com

 : 732-321-6915