




OUR OPTICS YOUR VISION




Microscope Objectives Plan Achromatic DW-PD Model 2



Plan Achromatic DW-PD Model 2				
Magnification	40X	50X	60X	100X (Oil)
Numerical Aperture	0.65	0.95	0.85	1.25
Working Distance (mm)	0.65	0.21	0.18	0.16
Slide CT (mm)	0.17	0.17	0.17	0.17
Focus Length (mm)	4.29	1.6	2.8	1.92
BFL (mm)				
Conjugate Distance (mm)	195	195	195	195
Tube Length (mm)	160	160	160	160
Wavelength (nm)	400-700	400-700	400-700	400-700
Resolution (μm)	0.4	0.3	0.2	0.3
Physical Depth of Focus (μm)	0.7	0.46	0.3	0.4
Thread	W4/5' x 1/36'			

 : www.shanghai-optics.com

 : rfq@shanghai-optics.com

 : 732-321-6915