





Microscope Objectives ∞ Inverted

| ∞ Inverted | | | | |
|--------------------------------|-----------|------------|------------|------------|
| Magnification | 4X | 10X | 20X | 40X |
| Numerical Aperture | 0.1 | 0.25 | 0.4 | 0.6 |
| Working Distance (mm) | 17.3 | 10 | 5.1 | 2.1 |
| Slide CT (mm) | - | 1.2 | 1.2 | 1.2 |
| BFL (mm) | | | | |
| Conjugate Distance (mm) | ∞ | ∞ | ∞ | ∞ |
| Wavelength (nm) | | | | |
| Resolution (μm) | | | | |
| Spring | | | | √ |
| Phase contrast | | √ | √ | √ |

 : www.shanghai-optics.com

 : rfq@shanghai-optics.com

 : 732-321-6915